


circaid® juxtalite™ lower leg



circumference in centimeters								
size	S	M	M – X	L	L – X	XL	XL – X	XXL
cC	25 – 36	33 – 44	44 – 54	41 – 51	51 – 61	37 – 48	48 – 61	53 – 64
cB	19 – 22	22 – 27	22 – 27	25 – 32	25 – 32	30 – 37	30 – 37	35 – 42
actual garment length in centimeters								
size	short				long			
ℓ B – D	28				33			
sizing measurement	short				long			
ℓ floor–knee crease	<44				>44			


All garments contain the circaid® built-in pressure system™ (BPS™) guide card.

Garment lengths may differ from body measurements.


circaid® capri

SIZE	Circumference in Inches		Circumference in Centimeters	
	WAIST (T)	THIGH (F)	WAIST (T)	THIGH (F)
A	25 - 31	18 - 24	63 - 79	46 - 62
B	26 - 33	19 - 26	67 - 85	49 - 65
C	28 - 36	20 - 27	72 - 91	52 - 70
D	30 - 38	22 - 29	76 - 96	55 - 73
E	32 - 40	24 - 32	81 - 102	60 - 81
F	33 - 43	25 - 33	85 - 108	63 - 84
G	35 - 44	26 - 35	90 - 113	67 - 90
H	37 - 47	27 - 37	94 - 119	70 - 93
I	39 - 49	28 - 38	99 - 125	72 - 96
J	41 - 51	30 - 40	103 - 130	77 - 102

When taking measurements, measure at the largest circumference point. When selecting a size, choose the smallest size match for the most snug fit. Choose the largest size match for a more relaxed, easy-to-wear to fit.



circaid® power added
compression band™ (pac band)
circaid® compression anklet



	circumferences in cm.	
	standard	large
cB	19 - 27	25 - 50
cH	19 - 36	22 - 41